

GC-QMS Application: Integrated Qualitative Analysis of Additive Components in Adhesive by using a Py-GC-QMS and msFineAnalysis iQ

Product used: Mass Spectrometer (MS)

Introduction

The gas chromatograph-quadrupole mass spectrometer (GC-QMS) is widely used as a qualitative / quantitative analyzer for volatile compounds. Qualitative analysis by GC-QMS is generally performed by searching the library database (DB) using the measurement data of the Electron Ionization (EI). However, if qualitative analysis is performed using only the similarity with the library spectrum as an index, several significant candidates may be obtained depending on the compound, or erroneous candidates may be selected as the identification result. In such cases, confirmation of molecular ions by the soft ionization (SI) such as the photoionization (PI) is effective.

In 2021, we released msFineAnalysis iQ, an integrated qualitative analysis software that automatically combines the analysis results of EI and SI methods measured by GC-QMS. The details of this software are introduced in MSTips No. 347 and 348. In this MSTips, we will introduce an analysis example of GC-MS measurement results for additive components in vinyl acetate adhesive using msFineAnalysis iQ.

Experimental

A commercially-available vinyl acetate adhesive was used as a test sample in this study. A GC-QMS (JMS-Q1600GC UltraQuad™ SQ-Zeta, JEOL Ltd.) was used for the measurement. We performed Py-GC-QMS measurements using both EI and photoionization (PI) modes with a combination EI/PI ion source. The qualitative data processing was performed with msFineAnalysis iQ (JEOL Ltd.). Detailed measurement conditions are shown in Table 1.



JMS-Q1600GC UltraQuad™ SQ-Zeta

Table 1 Measurement condition

Py (EGA/PY-3030D, Frontier Labs)		MS	
Sample amount	EI: 0.4 mg, PI: 1.0 mg	Ion Source Temp.	250°C
Pyrolysis Temp.	600°C	Interface Temp.	250°C
GC		Ion Source	EI/PI combination ion source
Column	ZB-5MS (Phenomenex) 30 m×0.25 mm I.D., df=0.25 μm	Ionization	EI+ (70 eV, 50 μA), PI+ (8~10 eV)
Injector Temp.	300°C	Acquisition Mode	Scan (m/z 35 - 600)
Oven Temp.	40°C (2 min) → 10°C/min → 320°C (5min)		
Injection Mode	Split 100:1		
Carrier Gas	He, 1.0 mL/min (Constant Flow)		

Results and Discussion

Figure 1 shows TICC of Py-GC/EI and Py-GC/PI measurement results. Strong peaks derived from Benzene, Acetic acid, Toluene, Indene, Naphthalene etc. were observed. In addition, components presumed to be additive components were also detected (components A, B and C). Detailed analysis results for components A, B and C are shown in the next section.




Figure 1 Py-GC/EI and Py-GC/PI total ion current chromatograms

Figure 2 shows mass spectra of components A, B and C (in EI mass spectrum, black line: measured spectrum, red line: library spectrum). For these components, molecular ions were detected in the mass spectra of both the EI and PI data (IM indicate molecular ion in Figure 2). However, the relative intensities of molecular ions were higher in the PI method. The integrated qualitative analysis result lists of the top 3 candidates each components is shown in Table 2. Component A was estimated as "Ethanol, 2-phenoxy-" because the similarity score calculated with the library DB was 922. The retention index value was Δ RI 2iu, which was a good value. This is a representative component used as a preservative in cosmetics and pharmaceuticals. In resins, it is used as a plasticizer and a film formation assistant. From the above results, it was found that additives were contained in the vinyl acetate adhesive. Component B was estimated as "Ethanol, 2-(2-phenoxyethoxy)-" because the similarity score calculated with the library DB was 935. The structural formula in Figure 2 shows that component B has the structure of the ethylene glycol moiety of component A (OCH₂CH₂OH) + ethylene oxide (OCH₂CH₂). And component C was estimated as "Ethanol, 2-[2-(2-phenoxyethoxy)ethoxy]-" because the similarity score calculated with the library DB was 917. The structural formula of component C was component B + ethylene oxide (OCH₂CH₂). These results indicate that components B and C are similar in structure to the additive component A. The similarity in structure suggests that the properties of the compounds may also be similar.




Figure 2 Mass spectra of component A,B and C

Table 2 Integrated qualitative analysis result of component A,B and C

[A]	#	Library Name	CAS#	Lib.	Similarity	Reverse Similarity	Lib. RI [iu]	Δ RI [iu]	Formula	DBE	El Base Peak (Lib.)	MW	Molecular Weight Check	Adduct /Loss	Isotope Matching
	L01	Ethanol, 2-phenoxy-	122-99-6	mainlib	922	925	1225	2	C8 H10 O2	4.0	94	138	✓	none	0.99
	L02	Carbamic acid, butylmethyl-, phenyl este	54644-61-0	mainlib	731	864	1126-1888	0	C12 H17 N O2	5.0	94	207	-	-	-
	L03	Carbonic acid, phenyl propyl ester	13183-16-9	mainlib	718	820	954-1716	0	C10 H12 O3	5.0	94	180	-	-	-
	L04	Carbonic acid, ethyl phenyl ester	3878-46-4	mainlib	702	799	1233	10	C9 H10 O3	5.0	94	166	-	-	-
	L05	1-Phenoxy-2-chloropropane	53491-30-8	mainlib	701	790	828-1590	0	C9 H11 Cl O	4.0	94	170	-	-	-
[B]	#	Library Name	CAS#	Lib.	Similarity	Reverse Similarity	Lib. RI [iu]	Δ RI [iu]	Formula	DBE	El Base Peak (Lib.)	MW	Molecular Weight Check	Adduct /Loss	Isotope Matching
	L01	Ethanol, 2-(2-phenoxyethoxy)-	104-68-7	mainlib	935	935	1105-1867	0	C10 H14 O3	4.0	45	182	✓	none	0.99
	L24	Benzene, 4-methyl-1,2-dinitro-	610-39-9	mainlib	531	555	1585	73	C7 H6 N2 O4	6.0	30	182	✓	none	0.81
	L02	Benzene, 1,1'-[1,2-ethanediylbis(oxy-2,1-ethanediylloxy)]bis-	53129-28-5	mainlib	640	648	1975-2559	463	C18 H22 O4	8.0	77	302	-	-	-
	L03	2-Phenoxyethanol, isopropyl ether	-	mainlib	625	686	987-1571	0	C11 H16 O2	4.0	94	180	-	-	-
	L04	2-Phenoxyethanol, n-propyl ether	-	mainlib	620	679	1051-1635	0	C11 H16 O2	4.0	45	180	-	-	-
[C]	#	Library Name	CAS#	Lib.	Similarity	Reverse Similarity	Lib. RI [iu]	Δ RI [iu]	Formula	DBE	El Base Peak (Lib.)	MW	Molecular Weight Check	Adduct /Loss	Isotope Matching
	L01	Ethanol, 2-[2-(2-phenoxyethoxy)ethoxy]-	7204-16-2	mainlib	917	936	1380-2142	0	C12 H18 O4	4.0	45	226	✓	none	0.95
	L02	Benzene, 1,1'-[oxybis(2,1-ethanediylloxy-2,1-ethanediylloxy)]bis-	20768-77-8	mainlib	706	724	2278-2862	480	C20 H26 O5	8.0	77	346	-	-	-
	L03	Ethanol, 2-(2-phenoxyethoxy)-	104-68-7	mainlib	646	686	1105-1867	0	C10 H14 O3	4.0	45	182	-	-	-
	L04	Benzene, 1,1'-[1,2-ethanediylbis(oxy-2,1-ethanediylloxy)]bis-	53129-28-5	mainlib	583	621	1975-2559	177	C18 H22 O4	8.0	77	302	-	-	-
	L05	2-Nitrophenethyl alcohol, methyl ether	-	mainlib	578	637	1083-1845	0	C9 H11 N O3	5.0	45	181	-	-	-

Conclusions

In this MSTips, we introduced analysis results of additive components in vinyl acetate adhesive using msFineAnalysis iQ. msFineAnalysis iQ automatically performs an integrated analysis that combines library DB search results using GC/EI data and confirmation of molecular ions using GC/SI data. Therefore, it was possible to easily analyze the additives in the polymeric material. This software is expected to improve the qualitative accuracy and efficiency of qualitative analysis using GC-QMS.

Copyright © 2023 JEOL Ltd.

Certain products in this brochure are controlled under the "Foreign Exchange and Foreign Trade Law" of Japan in compliance with international security export control. JEOL Ltd. must provide the Japanese Government with "End-user's Statement of Assurance" and "End-use Certificate" in order to obtain the export license needed for export from Japan. If the product to be exported is in this category, the end user will be asked to fill in these certificate forms.

